

Product Datasheet

PIGW peptide (orb218764)

Description	The peptide is used to block anti-PIGW antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PIGW
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215408.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7TSN4 , Q7Z7B1 , Q8C398
Entrez	70325 , 378774 , 284098
Expiration Date	12 months from date of receipt.