

Product Datasheet

GPR153 peptide (orb218778)

Description	The peptide is used to block anti-GPR153 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GPR153
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215422.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6NV75, Q8K0Z9
Entrez	100129, 387509
Expiration Date	12 months from date of receipt.