

Product Datasheet

SLC6A17 peptide (orb218780)

Description	The peptide is used to block anti-SLC6A17 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SLC6A17
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215424.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31662 , Q8BJI1 , Q9H1V8
Entrez	229706 , 388662 , 613226
Expiration Date	12 months from date of receipt.