

Product Datasheet

CD52 peptide (orb218808)

Description	The peptide is used to block anti-CD52 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CD52
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215452
Purity	> 85%
Source	Synthetic
Uniprot ID	Q64389, Q63064, P31358
Entrez	23833, 1043, 117054
Expiration Date	12 months from date of receipt.