

Product Datasheet

CLCN1 peptide (orb218810)

Description	The peptide is used to block anti-CLCN1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CLCN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215454.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35523, Q64347, P35524
Entrez	12723, 1180, 25688
Expiration Date	12 months from date of receipt.