

Product Datasheet

COX8A peptide (orb218811)

Description	The peptide is used to block anti-COX8A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	COX8A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215455.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10176, P80433, Q64445
Entrez	171335, 12868, 1351
Expiration Date	12 months from date of receipt.