

Product Datasheet

DYNC1I2 peptide (orb218816)

Description	The peptide is used to block anti-DYNC1I2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	DYNC1I2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215460.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13409, Q62871, O88487
Entrez	1781, 13427
Expiration Date	12 months from date of receipt.