

Product Datasheet

Goat Mouse Ig Lambda Light Chain (free and bound) antibody (orb21882)

Description

Precipitating polyclonal goat

Species/Host

Goat

Conjugation

Unconjugated

Preservatives

Delipidated, heat inactivated lyophilized whole serum. No preservative added, as it may interfere with the antibody activity

Form/Appearance

Delipidated, heat inactivated lyophilized whole serum. No preservative added, as it may interfere with the antibody activity

Storage

Storage: The lyophilized product is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

Note

For research use only

Clonality

Polyclonal

Source

A pool of purified Bence Jones lambda proteins isolated from mouse urine. Freund's complete adjuvant is used in the first step of the immunization procedure.

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this