

Product Datasheet

FCAR peptide (orb218824)

Description	The peptide is used to block anti-CD89 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FCAR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215468.
Purity	> 85%
Source	Synthetic
Uniprot ID	P24071
Entrez	2204
Expiration Date	12 months from date of receipt.