

Product Datasheet

GJB5 peptide (orb218825)

Description	The peptide is used to block anti-Connexin 31.1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GJB5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215469.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02739, O95377, P28232
Entrez	14622, 29586, 2709
Expiration Date	12 months from date of receipt.