

Product Datasheet

CXCR1 peptide (orb218835)

Description	The peptide is used to block anti-CD181 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CXCR1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215479.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70612, P25024, Q810W6
Entrez	54258, 227288, 3577
Expiration Date	12 months from date of receipt.