

Product Datasheet

MYH9 peptide (orb218851)

Description	The peptide is used to block anti-Myosin 9 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	MYH9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215495.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62812, Q8VDD5, P35579
Entrez	4627, 17886, 25745
Expiration Date	12 months from date of receipt.