

Product Datasheet

E-selectin peptide (orb218873)

Description	E-selectin peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SELE
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215517.
Purity	> 85%
Source	Synthetic
Uniprot ID	P98105, Q00690, P16581
Entrez	25544, 6401, 20339
Expiration Date	12 months from date of receipt.