

## Product Datasheet

### SLC9A5 peptide (orb218876)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | SLC9A5 peptide                                           |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rat                      |
| <b>Tested Applications</b> | Blocking                                                 |
| <b>Target</b>              | SLC9A5                                                   |
| <b>Form/Appearance</b>     | Lyophilized powder                                       |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.             |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb215520. |
| <b>Purity</b>              | > 85%                                                    |
| <b>Source</b>              | Synthetic                                                |
| <b>Uniprot ID</b>          | <a href="#">Q9Z0X2</a> , <a href="#">Q14940</a>          |
| <b>Entrez</b>              | <a href="#">192215</a> , <a href="#">6553</a>            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                          |