

Product Datasheet

SNAIL2 peptide (orb218877)

Description	The peptide is used to block anti-SLUG antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	SNAIL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215521.
Purity	> 85%
Source	Synthetic
Uniprot ID	O43623, P97469, O08954
Entrez	20583, 25554, 6591
Expiration Date	12 months from date of receipt.