

Product Datasheet

SSTR2 peptide (orb218880)

Description	The peptide is used to block anti-SSTR2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	SSTR2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215524.
Purity	> 85%
Source	Synthetic
Uniprot ID	P30680, P30874, P30875
Entrez	6752, 20606, 54305
Expiration Date	12 months from date of receipt.