

Product Datasheet

MED22 peptide (orb218882)

Description	The peptide is used to block anti-MED22 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MED22
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215526.
Purity	> 85%
Source	Synthetic
Uniprot ID	A0JPN6 , Q15528 , Q62276
Entrez	20933 , 499762 , 100911917 , 6837
Expiration Date	12 months from date of receipt.