

Product Datasheet

TERT peptide (orb218885)

Description	The peptide is used to block anti-Telomerase antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TERT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215529.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14746, O70372, Q673L6
Entrez	301965, 7015, 21752
Expiration Date	12 months from date of receipt.