

Product Datasheet

AXIN1 peptide (orb218888)

Description	The peptide is used to block anti-Axin-1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	AXIN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215532.
Purity	> 85%
Source	Synthetic
Uniprot ID	035625, 070239, 015169
Entrez	79257, 12005, 8312
Expiration Date	12 months from date of receipt.