

Product Datasheet

GEMIN2 peptide (orb218889)

Description	The peptide is used to block anti-Gemin2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GEMIN2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215533.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14893, Q9CQQ4, Q9QZP1
Entrez	8487, 84404, 66603
Expiration Date	12 months from date of receipt.