

Product Datasheet

SPAG9 peptide (orb218894)

Description	The peptide is used to block anti-JIP4 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SPAG9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215538.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q58A65, O60271
Entrez	9043, 70834
Expiration Date	12 months from date of receipt.