

Product Datasheet

LATS1 peptide (orb218895)

Description	The peptide is used to block anti-LATS1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	LATS1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215539.
Purity	> 85%
Source	Synthetic
Uniprot ID	O95835, Q8BYR2
Entrez	16798, 9113
Expiration Date	12 months from date of receipt.