

Product Datasheet

RHOQ peptide (orb218913)

Description	The peptide is used to block anti-TC10 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RHOQ
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215557.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8R527 , P17081 , Q9JJL4
Entrez	23433 , 104215 , 85428
Expiration Date	12 months from date of receipt.