

Product Datasheet

KAT6B peptide (orb218914)

Description	The peptide is used to block anti-MYST4 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	KAT6B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215558.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BRB7 , Q8WYB5
Entrez	23522 , 54169
Expiration Date	12 months from date of receipt.