

Product Datasheet

Goat Human IgA (Fc specific) antibody (orb21892)

Description

IgG fraction of polyclonal goat

unlabeled primary or secondary...

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

IP

Applications
Storage

Storage: The lyophilized antiserum is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Dilutions may be prepared by adding phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the antiserum. Diluted antiserum should be stored at +4°C, not refrozen, and preferably used the same day

Note

For research use only

Clonality

Polyclonal

Source

Highly purified normal IgA and monoclonal IgA2 isolated from human serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

Hazard
Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

Expiration Date

12 months from date of receipt.