

Product Datasheet

POLR3K peptide (orb218926)

Description	The peptide is used to block anti-POLR3K antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	POLR3K
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215570.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CQZ7, Q9Y2Y1
Entrez	67005, 51728
Expiration Date	12 months from date of receipt.