

Product Datasheet

GNG12 peptide (orb218933)

Description	The peptide is used to block anti-GNG12 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	GNG12
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215577.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9DAS9, Q9UBI6
Entrez	55970, 14701
Expiration Date	12 months from date of receipt.