

Product Datasheet

MKL2 peptide (orb218938)

Description	MKL2 peptide
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	MKL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb215582.
Purity	> 85%
Source	Synthetic
Uniprot ID	P59759, Q9ULH7
Entrez	57496, 239719
Expiration Date	12 months from date of receipt.