

Product Datasheet

HLA ABC Antibody - FITC Conjugated (orb2189542)

Description	HLA ABC Antibody - FITC Conjugated
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Dendritic cells from synovial fluid.
Preservatives	Liquid. Phosphate buffered saline
Form/Appearance	Liquid. Phosphate buffered saline
Storage	-20°C
Note	For research use only
Application notes	Application Info: Flow Cytometry: Â Â Neat - 1/10
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P30443
NCBI	NP_001229687.1
Expiration Date	12 months from date of receipt.