

Product Datasheet

CD84 Antibody - FITC Conjugated (orb2189665)

Description	CD84 Antibody - FITC Conjugated
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Mouse pre-B cell line 300.19 transfected with CD84 cDNA.
Preservatives	Liquid. Phosphate buffered saline
Form/Appearance	Liquid. Phosphate buffered saline
Storage	-20°C
Note	For research use only
Application notes	Application Info: Flow Cytometry: Â Â Neat
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9UIB8
NCBI	NP_001171808.1
Dilution Range	Application #1: Staining of human peripheral blood monocytes with Mouse anti Human CD84 Application #2: Staining of human peripheral blood monocytes with Mouse anti Human CD84: RPE
Expiration Date	12 months from date of receipt.