

Product Datasheet

Aquaporin 3 Peptides (orb218996)

Description	Aquaporin 3 Peptides
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	AQP3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216020.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q92482, Q8R2N1, P47862
Entrez	360, 65133, 11828
Expiration Date	12 months from date of receipt.