

Product Datasheet

Collagen IV peptide (orb219039)

Description	Collagen IV peptide
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	COL4A3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216063.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q01955, Q9QZS0
Entrez	12828, 1285
Expiration Date	12 months from date of receipt.