

Product Datasheet

GJA4 peptide (orb219094)

Description	The peptide is used to block anti-GJA4 reactivity
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	GJA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216118.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35212, P28235, Q03190
Entrez	2701, 14612, 25655
Expiration Date	12 months from date of receipt.