

Product Datasheet

LIMK1; LIMK2 peptide (orb219119)

Description	LIMK1; LIMK2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	LIMK1; LIMK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216143.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P53667, P53668, P53669
Entrez	16885, 3984
Expiration Date	12 months from date of receipt.