

Product Datasheet

MYF6 peptide (orb219140)

Description	The peptide is used to block anti-MYF6 reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MYF6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216164.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15375, P19335, P23409
Entrez	17878, 25714, 4618
Expiration Date	12 months from date of receipt.