

Product Datasheet

MAPK1 peptide (orb219161)

Description	The peptide is used to block anti-MAPK1 reactivity
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	MAPK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb216185.Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P63085, P63086, P28482
Entrez	116590, 26413, 5594
Expiration Date	12 months from date of receipt.