

## Product Datasheet

### Protein A Agarose IP Reagent [Out of stock] (orb219211)

<b>Description</b>	Protein A Agarose IP Reagent
<b>Tested Applications</b>	IP
<b>Storage</b>	Store at 4°C for one year. Do not freeze.
<b>Note</b>	For research use only
<b>Application notes</b>	Protein A is provided as an agarose conjugate for use in immunoprecipitation only. The product is provided as 0.5 ml agarose in 2.0 ml PBS. Protein A-Agarose is pre-blocked with BSA to reduce non-specific immunoglobulin binding. Sufficient product is provided for 100 immunoprecipitation reactions, to be used at 20 ul resuspended volume per reaction.
<b>Expiration Date</b>	12 months from date of receipt.