

Product Datasheet

Protein G Agarose IP Reagent [Out of stock] (orb219212)

Description	Protein G Agarose IP Reagent
Tested Applications	IP
Storage	Store at 4°C for one year. Do not freeze.
Note	For research use only
Application notes	Protein G is provided as an agarose conjugate for use in immunoprecipitation only. The product is provided as 0.5 ml agarose in 2.0 ml PBS. Protein G-Agarose is pre-blocked with BSA to reduce non-specific immunoglobulin binding. Sufficient product is provided for 100 immunoprecipitation reactions, to be used at 20 ul resuspended volume per reaction.
Expiration Date	12 months from date of receipt.