


Product Datasheet

Human Beta Amyloid ELISA Kit (orb219233)

Description	Sandwich ELISA for quantitative detection of Human Beta Amyloid concentrations.
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	APP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P05067
Sensitivity	40 pg/ml
Entrez	351
Expiration Date	Please enquire.



Representative
Standard
Curve