

## Product Datasheet

### Albumin antibody (orb21924)

**Description**

Goat polyclonal antibody to Albumin and is

albumin as te...

**Species/Host**

Goat

**Conjugation**

Unconjugated

**Preservatives**

Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those in pooled goat serum. No foreign proteins added.

**Form/Appearance**

Delipidated, heat inactivated, lyophilized, stable whole antiserum. No preservative added. Total protein and IgG concentrations in the antiserum are comparable to those in pooled goat serum. No foreign proteins added.

**Storage**

Storage: The lyophilized antiserum is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Dilutions may be prepared by adding phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the antiserum. Diluted antiserum should be stored at +4°C, not refrozen, and preferably used the same day

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Highly purified albumin isolated from cat serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

**Expiration Date**

12 months from date of receipt.