


Product Datasheet

Human Carbonic Anhydrase 9 ELISA Kit (orb219240)

Description	Sandwich ELISA for quantitative detection of Human Carbonic Anhydrase 9 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CA9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	Q16790
Sensitivity	10 pg/ml
Entrez	768
Expiration Date	Please enquire.



Representative
Standard
Curve