


Product Datasheet

Human TNFR1 ELISA Kit (orb219271)

Description	Human TNFR1 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	TNFRSF1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P19438
Sensitivity	1 pg/ml
Entrez	7132
Expiration Date	Please enquire.



Representative
Standard
Curve