

Product Datasheet

Human CD127 ELISA Kit (orb219274)

Description Sandwich ELISA for quantitative detection of Human CD127 concentrations.

Reactivity Human

Range 312 pg/ml - 20000 pg/ml

Tested Applications ELISA

Target IL7R

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Cell Culture Supernatants, Serum
Experiment Notes: Colorimetric


Sample Types Cell Culture Supernates, Serum

Uniprot ID **P16871**

Sensitivity 10 pg/ml

Entrez **3575**

Expiration Date Please enquire.



Representative
Standard
Curve