


Product Datasheet

Human CD257 ELISA Kit (orb219298)

Description	Sandwich ELISA for quantitative detection of Human CD257 concentrations.
Reactivity	Human
Range	39 pg/ml - 2500 pg/ml
Tested Applications	ELISA
Target	TNFSF13B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9Y275
Sensitivity	20 pg/ml
Entrez	10673
Expiration Date	Please enquire.



Representative
Standard
Curve