

## Product Datasheet

### Lysozyme antibody (orb21930)

<b>Description</b>	Rabbit polyclonal antibody to Lysozyme.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Immunogen</b>	Highly purified recombinant human lysozyme was used as immunogen.
<b>Storage</b>	Storage: The lyophilized IgG fraction is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Lysozyme, a bacteriolytic enzyme, has a molecular weight of about 15,000. Its concentration in plasma is low; in milk the concentration is reported to be between 5 and 220 mg/ml. Also tears contain a considerable amount of lysozyme. Highly purified recombinant human lysozyme was used as immunogen. Freund's complete adjuvant is used in the first step of the immunization procedure.
<b>Uniprot ID</b>	<b>P61626</b>
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.
<b>Expiration Date</b>	12 months from date of receipt.