


Product Datasheet

Human RANKL ELISA Kit (orb219301)

Description	Human RANKL ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	TNFRSF11A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	Q9Y6Q6
Sensitivity	2 pg/ml
Entrez	8792
Expiration Date	Please enquire.



Representative
Standard
Curve