


Product Datasheet

Human CD27 ELISA Kit (orb219305)

Description	Sandwich ELISA for quantitative detection of Human CD27 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	CD27
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, SerumExperiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	P26842
Sensitivity	10 pg/ml
Entrez	939
Expiration Date	Please enquire.



Representative
Standard
Curve