


Product Datasheet

Human CD62P ELISA Kit (orb219323)

Description	Human CD62P ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	SELP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P16109
Sensitivity	30 pg/ml
Entrez	6403
Expiration Date	Please enquire.



Representative
Standard
Curve