


Product Datasheet

Human CD66e ELISA Kit (orb219324)

Description	Sandwich ELISA for quantitative detection of Human CD66e concentrations.
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	CEACAM5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P06731
Sensitivity	40 pg/ml
Entrez	1048
Expiration Date	Please enquire.



Representative
Standard
Curve