


Product Datasheet

Human CD87 ELISA Kit (orb219326)

Description	Sandwich ELISA for quantitative detection of Human CD87 concentrations.
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	PLAUR
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Cell Culture Supernatants, Serum, Plasma, Urine Experiment Notes: Colorimetric
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	Q03405
Sensitivity	4 pg/ml
Entrez	5329
Expiration Date	Please enquire.



Representative
Standard
Curve